Guardian™ 4

Sensor User Guide

Korisnički priručnik za senzor

Korisnički priručnik za senzor

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Introduction

The Guardian™ 4 sensor is part of the Continuous Glucose Monitoring (CGM) system. The sensor converts small amounts of glucose from the interstitial fluid under the skin into an electronic signal. The system then uses these signals to provide sensor glucose values.

Sensor assembly

Indications for use

The Guardian 4 sensor (MMT-7040Q) is intended for use with the Guardian™ 4 transmitter (MMT-7841Q) to monitor glucose levels in persons with diabetes. The sensor is designed to replace fingerstick blood glucose readings for diabetes treatment decisions.

The sensor is intended for insertion into persons ages 7 years and older. The sensor is intended for insertion into the back of the upper arm or the upper buttocks in persons ages 7 through 17 years. The sensor is intended for insertion into the back of the upper arm or the abdomen in persons ages 18 years and older.

Contraindications

No contraindications are associated with Guardian 4 sensor use. For contraindications related to CGM, see the Guardian™ App User Guide.

Clinical benefits

The Guardian 4 sensor is a component of the CGM system that provides sensor glucose values. See the Guardian App User Guide for the clinical benefits of systems that use the Guardian 4 sensor.

User safety

Warnings

Read this entire user guide before attempting to insert the Guardian 4 sensor. The one-pressserter (MMT-7512) is the onlyserter approved for use with the sensor. Failure to follow directions, or the use of a different insertion device, may result in improper insertion, pain, or injury.

Do not attempt to connect a transmitter or recorder that is not compatible with the sensor. The sensor is designed to work with approved transmitters only. Connecting the sensor to a transmitter or recorder that is not approved for use with the sensor may damage the components. Refer to the Guardian App User Guide for a list of compatible products.

Do not use continuous glucose monitoring if hydroxyurea, also known as hydroxycarbamide, is taken. Hydroxyurea is used to

treat certain diseases, such as cancer and sickle cell anemia. Hydroxyurea use results in higher sensor glucose readings compared to blood glucose readings. Taking hydroxyurea while using continuous glucose monitoring can result in substantially higher sensor glucose readings in reports than actual blood glucose readings.

Always check the label of any medication being taken to confirm if hydroxyurea or hydroxycarbamide is an active ingredient. If hydroxyurea is taken, consult a healthcare professional. Use additional blood glucose meter readings to verify glucose levels.

Taking medications that contain acetaminophen or paracetamol, including, but not limited to fever reducers and cold medicine, while wearing the sensor, may falsely raise sensor glucose readings. The level of inaccuracy depends on the amount of acetaminophen or paracetamol active in the body and may be different for each person. Always check the label of any medications to confirm whether acetaminophen or paracetamol is an active ingredient.

Do not expose the sensor to MRI equipment, diathermy devices, or other devices that generate strong magnetic fields. The performance of the sensor has not been evaluated under those conditions and may be unsafe. If the sensor is exposed to a strong magnetic field, discontinue use and contact a local Medtronic support representative for further assistance.

Always inspect the packaging for damage prior to use. Sensors are sterile and non-pyrogenic, unless the package has been opened or damaged. If the sensor packaging is open or damaged, discard the sensor directly into a sharps container. Use of a non-sterile sensor may result in infection at the insertion site.

Do not allow children to put small parts in their mouth. This product may pose a choking hazard that can result in serious injury or death.

Healthcare professionals and caregivers:

- Always wear gloves to insert the sensor. A retractable needle is attached to the sensor. Minimal bleeding may occur.
- Cover the sensor with sterile gauze to remove the needle housing from the sensor.

Place the needle housing directly into a sharps container after sensor insertion to prevent accidental needlestick injury.

Watch for bleeding at the insertion site (under, around, or on top of the sensor).


If bleeding occurs, do the following:

1. Apply steady pressure, using sterile gauze or a clean cloth placed on top of the sensor, for up to three minutes. The use of unsterile gauze may cause site infection.
2. If bleeding stops, connect the transmitter (or recorder) to the sensor.

If bleeding does not stop, do not connect the transmitter to the sensor because blood may get into the transmitter connector, and may damage the device.

If bleeding continues, causes excessive pain or discomfort, or is significantly visible in the plastic base of the sensor, do the following:

1. Remove the sensor and continue to apply steady pressure until the bleeding stops. Discard the sensor in a sharps container.
2. Check the site for redness, bleeding, irritation, pain, tenderness, or inflammation. Treat based on instructions from a healthcare professional.
3. Insert a new sensor in a different location.

Plastic base

For questions or concerns related to sensor use, contact a local Medtronic support representative for assistance.

For medical questions or concerns, contact a healthcare professional.

If a serious incident related to the device occurs, immediately report the incident to a healthcare professional.

Healthcare professionals: report serious incidents to Medtronic and the appropriate authorities.

Precautions

Wash hands with soap and water before inserting the Guardian 4 sensor to help prevent site infection.

Do not insert the sensor through tape. Inserting the sensor through tape may cause improper sensor insertion and function.

Only use alcohol to prepare the insertion site. Using alcohol to prepare the insertion site makes sure that residue is not left on the skin.

Rotate the sensor insertion site so that sites do not become overused.

Do not clean, resterilize, or try to extract the needle from the needle housing. An accidental needlestick or puncture may occur.

Do not reuse sensors. Reuse of a sensor may cause damage to the sensor surface and lead to inaccurate glucose values, site irritation, or infection.

Risks and side effects

Keep the Guardian 4 sensor out of reach of children. The sensor may pose a choking hazard that can result in serious injury or death.

Other risks related to sensor use include:

- Skin irritation or other reactions
- Bruising
- Discomfort

- Redness
- Bleeding
- Pain
- Rash
- Infection
- Raised bump
- Appearance of a small “freckle-like” dot where needle was inserted
- Allergic reaction
- Fainting secondary to anxiety or fear of needle insertion
- Soreness or tenderness
- Swelling at insertion site
- Sensor fracture, breakage or damage
- Minimal blood splatter associated with sensor needle removal
- Residual redness associated with adhesive or tapes or both
- Scarring

Hazardous substances

None.

Allergens

None known.

Reagents






The Guardian 4 sensor contains two biological reagents: glucose oxidase, and human serum albumin (HSA). Glucose oxidase is derived from *Aspergillus niger* and manufactured to meet industry requirements for the extraction and purification of enzymes for use in diagnostic, immunodiagnostic, and biotechnical applications.

The HSA used on the sensor consists of purified and dried albumin fraction V derived from pasteurized human serum, which is cross-linked via glutaraldehyde. Approximately 3 µg of glucose oxidase and approximately 10 µg of HSA are used to manufacture each sensor. HSA is approved for IV infusion in humans at quantities much larger than in the sensor.

Remove the sensor

To change the Guardian 4 sensor, disconnect the transmitter from the sensor as described in the Guardian transmitter user guide. Gently pull the sensor from the body to remove it. Discard the sensor in a sharps container.

Components

	<p>One-pressserter A. bump on both buttons B. thumbprint marking</p>
	<p>Glucose sensor assembly A. pedestal B. needle housing C. sensor D. clear liner</p>
	<p>Sensor base A. sensor connector B. sensor snaps</p>
	<p>Transmitter</p>
	<p>First piece of oval tape A. liner 1 B. liner 2</p>
	<p>Second piece of oval tape A. liner 1 B. liner 2</p>

Where to insert the sensor

Choose an insertion site for the applicable age group. Target the shaded areas, and make sure that the insertion site has a sufficient amount of fat.

Ages 7-17 years

Upper buttocks

Back of the upper arm


Insertion on abdomen for ages 7-17 years has not been evaluated for accuracy.

Note: Insertion onto the upper buttocks should target the top third of the buttocks area. Assistance from another person may be needed for sensor insertion into the back of the upper arm or the upper buttock. If assistance is not needed, a mirror may be helpful for self insertion.

Ages 18 years and older

Back of the upper arm

Abdomen

Insertion on upper buttocks for ages 18 years and older has not been evaluated for accuracy.

CAUTION: Avoid the 5.0 cm (2 inch) area around the navel to help make a comfortable insertion site and to help with sensor adhesion.


For best sensor glucose performance, and to prevent accidental sensor removal:

- Do not insert the sensor into muscle, tough skin, or scar tissue.
- Avoid areas that are constrained by clothing or accessories.


- Avoid areas that are located under a belt or waistband.
- Avoid areas subjected to vigorous movement during exercise.

Inserting the sensor


WARNING: Always wear gloves when inserting the sensor into another person to avoid contact with patient blood. Minimal bleeding may occur. Contact with patient blood may cause infection.


1. Wash hands thoroughly with soap and water.


2. Choose an insertion site that has a sufficient amount of fat.

3. Clean the insertion site with alcohol. Let the area air dry.

4. Open the sensor package.

5. Hold the pedestal and remove the glucose sensor assembly from the package. Place the pedestal on a clean, flat surface such as a table.

tucked tab


Correct


Incorrect

6. Confirm that the adhesive tab of the sensor is tucked under the sensor connector and sensor snaps.



Correct


Incorrect

7. Using either hand, place a thumb on the thumbprint marking to hold the serter. Fingers must not touch the serter buttons.


thumb on thumbprint marking

8. Push the serter down onto the pedestal until the base of the serter sits flat on the table and there is a click.


9. With either hand, place two fingers on the base of the pedestal. With the other hand, grip the serter and pull the serter upwards.


Note: The arrow on the side of the serter aligns with the needle inside the serter.

arrow


WARNING: Never point a loaded serter toward any body part where insertion is not desired. An accidental button-push may cause the needle to inject the sensor in an undesired location, causing minor injury.


10. Place the serter on top of the prepared insertion site.

11. Press and release both serter buttons at the same time. Continue to hold the serter on top of the insertion site for five seconds or more to let the adhesive stick to the skin.


12. Lift the serter away from the insertion site. Fingers must not press the buttons while lifting the serter.

Sensor base

- A. sensor snaps
- B. sensor connector
- C. adhesive tab
- D. adhesive liner
- E. adhesive pad


If the sensor is inserted without assistance, complete step 13a. If a healthcare professional or caregiver assisted with sensor insertion, complete step 13b.

Patient:

13. a. Hold the sensor base against the skin at the sensor connector and at the opposite end of the sensor base. Hold the needle housing at the top and pull away from the sensor.


OR

Healthcare professional or caregiver:

13. b. Wrap sterile gauze around the sensor. Hold the sensor base against the skin at the sensor connector and the opposite end of the sensor base. Hold the needle housing at the top and pull away from the sensor.


WARNING: Always watch for bleeding at the insertion site. If bleeding occurs under, around, or on top of the sensor, apply steady pressure using sterile gauze or a clean cloth placed on top of the sensor for up to three minutes. The use of unsterile gauze may cause an infection. If bleeding does not stop, remove the sensor and apply steady pressure until the bleeding stops.

Note: After insertion, use of adhesive products, such as Skin Tac™, in addition to the oval tape is optional. If optional adhesive products are used, apply them to the skin under the adhesive pad prior to removing the liner. Adhesive products may also be applied to the adhesive pad or the skin around the sensor base. Allow the product to dry before continuing.

14. Remove the adhesive liner from under the adhesive pad. Pull the liner away from the sensor, staying close to the skin. Do not pull on the sensor when you remove the liner.

Note: Do not remove the adhesive liner from the rectangular adhesive tab. This tab will be used to secure the transmitter in a later step.

Note: If the sensor base moves, hold the sensor base down.

15. Firmly press the adhesive pad against the insertion site to confirm that the sensor base remains on the skin.


16. Untuck the adhesive tape from under the sensor connector.


17. Straighten the sensor adhesive tab so that it lies flat against the skin.

Applying Oval Tape


1. Remove the liner marked 1.

2. Apply the tape as shown and press down firmly.


wide part of tape covers half of sensor base

3


3. Remove liners marked 2 from each side.

4

4. Smooth the tape.


5

5. Connect the transmitter to the sensor.

Note: Wait for the green light on the transmitter to flash. If the green light does not flash, refer to the Troubleshooting section of the Guardian 4 transmitter user guide.

6

6. Cover the transmitter with the adhesive tab.


Note: Do not pull the tab too tightly.

7

7. To apply a second tape, remove liner marked 1.


8

8. Apply the second tape in the opposite direction to the first tape and place it on the transmitter. Press down firmly.


wide part of tape covers end of transmitter and skin

9

9. Remove liners marked 2 from each side.

10

10. Smooth the tape.

Note: For details on how to enter sensor settings into a compatible display device, refer to the Guardian App user guide.

Maintenance

Cleaning

The Guardian 4 sensor is a single-use, disposable device. No cleaning or maintenance is required.

Storage

CAUTION: Do not freeze the sensor, or store it in direct sunlight, extreme temperatures, or humidity. These conditions may damage the sensor.

Only store sensors at room temperature between 2 °C to 27 °C (36 °F to 80 °F).

Discard sensor after the “Use-by date” indicated on the label, if the package is damaged, or if the seal is broken.

Disposal

Dispose of the Guardian 4 sensor into a sharps container.


Technical specifications

Approximate dimensions
3.8 x 6.7 x 5.2 centimeters (1.5 x 2.6 x 2.0 inches)
Approximate weight
2.80 grams (.09 ounces)

Sensor life of use

The Guardian 4 sensor can be used one time and has a maximum life of 170 hours (seven days). The 170-hour life span of the sensor begins when the sensor is connected to the transmitter.

Icon table	
	Use-by date
	Medical device
	Importer
	Do not re-use
	Caution: consult instructions for use for important warnings or precautions not found on the label.
(1x)	One sensor per container/package
(5x)	Five sensors per container/package
(2x)	Two tapes per package
(10x)	Ten tapes per package
	Consult instructions for use
	Catalogue number
	Batch code

	Sterilized using irradiation
	Do not use if package is damaged
	Single sterile barrier system
	Storage temperature limit
	Open here
	Manufacturer
	Date of manufacture
	Do not resterilize
	Fragile, handle with care
	Keep dry
	Recyclable, contains recycled content
	Authorized representative in the European Community
	Conformité Européenne (European Conformity). This symbol means that the device fully complies with applicable European Union Acts.
	Magnetic Resonance (MR) unsafe
	Non-pyrogenic
	Requires prescription in the USA

©2021 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. Skin Tac™

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




Komponente


 <p>A diagram of a grey, dome-shaped device. Label A points to a blue rectangular area on the front. Label B points to a series of horizontal ridges on the side.</p>	<p>One-press server A. ispuščenje na oba dugmeta B. oznaka otiska prsta</p>
 <p>A diagram of a glucose sensor assembly. Label A points to the base. Label B points to the needle housing. Label C points to the sensor itself. Label D points to the protective cover.</p>	<p>Sastavni delovi senzora za glukozu A. postolje B. kućište igle C. senzor D. providna zaštita</p>
 <p>A diagram of the sensor base. Label A points to the sensor connector. Label B points to the sensor fasteners.</p>	<p>Osnova senzora A. konektor senzora B. pričvršćivači senzora</p>
 <p>A diagram of a small, white, oval-shaped transmitter.</p>	<p>Predajnik</p>
 <p>A diagram of an oval adhesive strip with a central hole. Label A points to the top adhesive layer. Label B points to the side adhesive layers.</p>	<p>Prvi komad ovalne trake A. obloga 1 B. obloga 2</p>
 <p>A diagram of a second oval adhesive strip, identical to the first one, with labels A and B.</p>	<p>Drugi komad ovalne trake A. obloga 1 B. obloga 2</p>

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Medtronic B.V.